● **TOYOBO ENZYMES** ●  
**(Diagnostic Reagent Grade)**

# ALKALINE PHOSPHATASE

*from Microorganism*

Orthophosphoric-monoester phosphohydrolase (alkaline optimum) (EC 3.1.3.1)

## PREPARATION and SPECIFICATION

Appearance	: Transparent liquid
Activity	: Grade II 30,000U/ml or more
Contaminants	: Adenosine deaminase $\leq 1.0 \times 10^{-4}\%$ Phosphodiesterase $\leq 3.0 \times 10^{-3}\%$

## PROPERTIES


Stability	: Stable at 4°C	
Molecular weight	: approx. 100,000	
Optimum pH	: 9.8	(Fig.1)
Optimum temperature	: $\geq 60^\circ\text{C}$	(Fig.2)
pH Stability	: pH 6.0–10.4 (25°C, 16hr)	(Fig.3)
Thermal stability	: below 60°C (pH 7.0, 60min)	(Fig.4)

## APPLICATIONS

This enzyme is useful for molecular biology.

## ASSAY

### Principle:

The appearance of p-Nitrophenol is measured at 405nm by spectrophotometry.

### Unit definition:

One unit causes the formation of one micromole of p-Nitrophenol per minute under the conditions described below.

### Method:

#### Reagents

- A. Diethanolamine buffer, pH10.25 : 5.0mM [Dilute 9.66ml of diethanolamine (MW=105.14) in 60ml of H<sub>2</sub>O, add 0.25ml of 0.1M MgCl<sub>2</sub> and, after adjusting the pH to 10.25 with 2N HCl, fill up to 100ml with H<sub>2</sub>O] (Prepare freshly)
- B. pNPP : 0.1M [Dissolve 371mg of p-nitrophenylphosphate disodium salt (MW=371.16) in 10ml buffer solution A] (Prepare freshly)
- C. Enzyme diluent : 0.1M Diethanolamine buffer pH 9.5 contg. 0.25mM MgCl<sub>2</sub>

#### Procedure

- Prepare the following reaction mixture in a cuvette (d=1.0cm) and equilibrate at 37°C for about 5 minutes.
 

Concentration in assay mixture	
Diethanolamine buffer	0.97 M
p-Nitrophenylphosphate	10 mM
MgCl <sub>2</sub>	0.25mM

2.6ml Buffer solution (A)

0.3ml Substrate solution (B)

- Add 0.1ml of the enzyme solution\* and mix by gentle inversion.
- Record the increase in optical density at 405nm against water for 3 to 5 minutes in a spectrophotometer thermostated at 37°C, and calculate the  $\Delta OD$  per minute from the initial linear portion of the curve ( $\Delta OD$  test).

At the same time, measure the blank rate ( $\Delta OD$  blank) by using the same method as the test except that the enzyme diluent is added instead of the enzyme solution.

\* Dilute the enzyme preparation to 0.1–0.3U/ml with ice-cold enzyme diluent (C), immediately before assay.

#### Calculation

Activity can be calculated by using the following formula :

$$\text{Volume activity (U/ml)} = \frac{\Delta OD/\text{min} (\Delta OD \text{ test} - \Delta OD \text{ blank}) \times V_t \times df}{18.5 \times 1.0 \times V_s} = \Delta OD/\text{min} \times 1.62 \times df$$


V<sub>t</sub> : Total volume (3.0ml)

V<sub>s</sub> : Sample volume (0.1ml)

18.5 : Millimolar extinction coefficient of p-Nitrophenol under the assay condition (cm<sup>2</sup>/micromole)


1.0 : Light path length (cm)

df : Dilution factor


**Fig.1. pH-Activity**

( in 1M Diethanolamine buffer, pH 8-10.5)


**Fig.3. pH-Stability**

( 25°C, 16hr-treatment with 0.1M buffer solution: pH 4-6, dimethylglutaric acid-NaOH; pH 6-8, K-phosphate; pH 8-9, Tris-HCl; pH 9-10, glycine-NaOH. Enzyme concentration: 10U/ml )

**Fig.2. Temperature activity**

( in 1M Diethanolamine buffer, pH 10.25)


**Fig.4. Thermal stability**

( 15min-treatment with 50mM K-phosphate buffer, pH 7.0. Enzyme concentration: 10U/ml )

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