

● **TOYOBO ENZYMES** ●
(Diagnostic Reagent Grade)

PEROXIDASE

from Horseradish

Donor:hydrogen-peroxidase oxidoreductase (EC 1. 11. 1. 7)



PREPARATION and SPECIFICATION

PEO-301/302 are isolated from horseradish root in our own method. PEO-131 is chromatographically separated from the other isoenzymes using SP Sephadex C-50 by the modified method of Paul et al.¹⁾ The peroxidase fractions having equal RZ values (ca.3.3) are combined to generate the preparation. This preparation has an RZ value of ca.3.3 and is electrophoretically homogeneous. On the other hand, Grade III is partially purified preparation.

Appearance	: Reddish-brown amorphous powder, lyophilized
Activity	: Grade I 250 Purpurogallin U/mg-solid or more (-131) (RZ ≥ 3.0, salt free) Grade III 110 Purpurogallin U/mg-solid or more (-301) (RZ ≥ 2.0, containing approx. 30% of stabilizers) Grade III 180 Purpurogallin U/mg-solid or more (-302) (RZ ≥ 2.0, salt free)
Contaminant	: Phosphatase ≤1.0×10 ⁻³ % (Grade III)

PROPERTIES

Stability	: Stable at -20°C for at least one year	(Fig.1,2,3)
Molecular weight	: approx. 40,000	
Structure	: Glycoprotein with one mole of protohaemin IX ²⁾	
Inhibitors	: Cyanide, sulfide, fluoride, azide ³⁾	
Optimum pH	: 6.0—7.0	(Fig.6)
Optimum temperature	: 45°C	(Fig.7)
pH Stability	: pH 5.0—10.0 (25°C, 20hr)	(Fig.8)
Thermal stability	: below 50°C (pH 6.0, 10min)	(Fig.9)
Effect of various chemicals	: (Table 1)	

APPLICATIONS^{4~11)}

This enzyme is useful for enzymatic determination of H₂O₂ in clinical analysis. Especially, the highly purified preparation (Grade I) is useful as a protein tracer in histo-and cyto-chemistry and as a valuable experimental tool in hodological neurography. Also, the enzyme preparation has been used as an enzyme label in enzyme immunoassay. Grade III (-302) is suitable for dry chemistry. On the other hand, the enzymes contribute for the reduction of phehol in waste water.

ASSAY

Principle:



The appearance of Purpurogallin is measured at 420nm by spectrophotometry.

Unit definition:

One purpurogallin unit causes the formation of one milligram of purpurogallin in 20 seconds under the conditions described below.

Method:

Reagents

- A. Pyrogallol solution : 5% (W/V)(Should be prepared fresh).
 B. H₂O₂ solution : 0.147M [Dilute 1.67ml of 30% (W/V) H₂O₂ to 100ml with H₂O] (Should be prepared fresh)
 C. Phosphate buffer, pH6.0 : 0.1M
 D. H₂SO₄ solution : 2.0N

Procedure

1. Prepare the following reaction mixture in a test tube (32 φ × 200mm) and equilibrate at 20°C for about 5 minutes.

14.0ml	H ₂ O	
2.0ml	Pyrogallol solution	(A)
1.0ml	H ₂ O ₂ solution	(B)
2.0ml	Phosphate buffer, pH6.0	(C)

Concentration in assay mixture	
Phosphate buffer	15 mM
Pyrogallol	40 mM
H ₂ O ₂	7.4mM

2. Add 1.0ml of the enzyme solution* and mix.
 3. After exactly 20 seconds at 20°C, add 1.0ml of 2.0 N H₂SO₄ solution (D) to stop the reaction.
 4. Extract the produced purpurogallin from the above stopped reaction mixture in five times with 15ml portions of ether and fill up the combined ether extracts to 100ml with fresh ether.
 5. Measure the optical density at 420nm against water (OD test).

At the same time, prepare the blank by first mixing the reaction mixture with 1.0ml of 2.0 N H₂SO₄ solution (D) after 20 a sec-incubation at 20°C, followed by the addition of the enzyme solution and extracting with ether by the same procedure as the test (OD blank).

- * Dissolve the enzyme preparation in ice-cold 0.1M phosphate buffer, pH 6.0 (C), dilute to 3.0–6.0 purpurogallin U/ml with the same buffer and store on ice.

Calculation

Activity** can be calculated by using the following formula :

$$\text{Volume activity (U/ml)} = \frac{\Delta \text{OD (OD test} - \text{OD blank)} \times \text{df}}{0.117 \times \text{Vs}} = \Delta \text{OD} \times 8.547 \times \text{df}$$

$$\text{Weight activity (U/mg)} = (\text{U/ml}) \times 1/\text{C}$$

Vs : Sample volume (1.0ml)

0.117: Optical density at 420 nm corresponding to 1mg% of Purpurogallin in ether.

df : Dilution factor

C : Enzyme concentration in dissolution (c mg/ml)

**One purpurogallin unit is equivalent to 13.5 international units determined with o-dianisidine at 25°C.

REFERENCES

- 1) K.G.Paul and T.Stigbrand; *Acta Chem.Scand.*, **24**, 3607 (1970).
- 2) L.M.Shannon et al.; *J.Biol.Chem.*, **241**, 2166 (1966).
- 3) E.Kay et al.; *J.Biol.Chem.*, **242**, 2470 (1967).
- 4) R.Lasek et al.; *Brain Res.*, **8**, 319 (1968).
- 5) W.M.Cowan et al.; *Brain Res*, **37**, 21 (1972).
- 6) J.H.La Vail et al.; *Brain Res.*, **58**, 470 (1973).
- 7) A.M.Graybiel and M.Devor; *Brain Res.*, **68**, 167 (1974).
- 8) A.H.Bunt et al.; *Brain Res.*, **102**, 152 (1976).
- 9) D.R.Colman et al.; *Brain Res.*, **102**, 156 (1976).
- 10) M.Dubois-Dalcq et al.; *J.Histochem.Cytochem.*, **25**, 1201 (1977).
- 11) M.Sato et al.; *Brain Res.*, **140**, 149 (1978).

Table 1. Effect of Various Chemicals on Peroxidase

[The enzyme dissolved in 0.1M phosphate buffer, pH 6.0 (50U/ml) was incubated with each chemical at 25°C for 1hr.]

Chemical	Concn.(mM)	*Residual activity(%)	Chemical	Concn.(mM)	*Residual activity(%)
None	—	100	MIA	2.0	99
Metal salt	2.0		NEM	2.0	97
MgCl ₂		102	IAA	2.0	99
CaCl ₂		102	Hydroxylamine	2.0	98
Ba(OAc) ₂		105	EDTA	5.0	95
FeCl ₃		98	o-Phenanthroline	2.0	98
CoCl ₂		97	α, α' -Dipyridyl	1.0	96
MnCl ₂		97	Borate	50	98
ZnCl ₂		99	NaF	2.0	98
CdCl ₂		99	NaN ₃	2.0	75
NiCl ₂		96	Triton X-100	0.10%	98
CuSO ₄		98	Brij 35	0.10%	80
Pb(OAc) ₂		96	Tween 20	0.10%	89
AgNO ₃		91	Span 20	0.10%	98
HgCl ₂		92	Na-cholate	0.10%	97
2-Mercaptoethanol	2.0	94	SDS	0.05%	98
PCMB	1.0	98	DAC	0.05%	102

*Residual activity was measured by 4AA-DEA method

4AA, 4-Aminoantipyrine; DEA, Diethylaniline

Ac, CH₃CO; PCMB, p-Chloromercuribenzoate; MIA, Monoiodoacetate; NEM, N-Ethylmaleimide; IAA, Iodoacetamide; EDTA, Ethylenediaminetetraacetate; SDS, Sodium dodecyl sulfate; DAC, Dimethylbenzylalkylammonium chloride.




Fig. 1. Stability (PEO-131) (Powder form)

[kept under dry conditions]




Fig. 2. Stability (PEO-301) (Powder form)

[kept under dry conditions]




Fig. 3. Stability (PEO-302) (Powder form)

[kept under dry conditions]




Fig. 4. Stability (Powder form)

[kept under dry conditions]




Fig. 5. Stability (Powder form)

[kept under dry conditions]




Fig. 6. pH-Activity

[20°C, 20sec-reaction in 0.1M buffer solution: pH4.0-6.0, acetate; pH6.0-8.0, phosphate]




Fig. 7. Temperature activity

[20sec-reaction in 0.1M phosphate buffer, pH6.0]




Fig. 8. pH-Stability

[25°C, 20hr-treatment with 50mM buffer solution: pH3.5-6.0, acetate; pH6.0-8.0, phosphate; pH9.0-11.0 borate]





Fig. 9. Thermal stability

[10min-treatment with 50mM phosphate buffer, pH6.0]

活性測定法 (Japanese)

1.原理



生成するPurpurogallinをエーテル抽出し、420nmの吸光度の変化で測定する。

2.定義

下記条件下で20秒間に1.0mgのPurpurogallinを生成する酵素量を1Purpurogallin単位(U)とする。

3.試薬

- A. 5%(W/V)ピロガロール水溶液(用時調製)
 - B. 0.147M H₂O₂水溶液 [30%(W/V)H₂O₂溶液 1.67mlを蒸留水で希釈して100mlとする] (用時調製)
 - C. 0.1Mリン酸緩衝液,pH6.0(反応混液及び酵素希釈用)
 - D. 2.0N H₂SO₄溶液
- 酵素溶液：酵素標品を予め氷冷した0.1Mリン酸緩衝液、pH6.0で溶解し、同緩衝液で3.0～6.0Purpurogallin U/mlに希釈して氷冷保存する。

4.手順

- ①試験管(32φ×200mm)に下記反応混液を調製し、20℃で約5分間予備加温する。

14.0ml	蒸留水	
2.0ml	ピロガロール水溶液	(A)
1.0ml	H ₂ O ₂ 水溶液	(B)
2.0ml	リン酸緩衝液	(C)
- ②酵素溶液1.0mlを加え、反応を開始する。
- ③20℃で正確に20秒間反応させた後、H₂SO₄溶液(D)1.0mlを加えて反応を停止させる。反応停止後の混液から生成したPurpurogallinをエーテル15mlで抽出する。この操作を5回繰り返し、抽出液を合わせ、更にエーテルを加えて全量を100mlにする。この液につき420nmにおける吸光度を測定する(OD test)。
- ④盲検は反応混液①を20℃で20秒間放置後、H₂SO₄溶液(D)1.0mlを加えて混和し、次いで酵素溶液1.0mlを加えて調製する。この液につき上記同様にエーテル抽出を行って吸光度を測定する(ODblank)。

5.計算式

$$\text{U/ml} = \frac{\Delta \text{OD} (\text{OD test} - \text{OD blank}) \times \text{希釈倍率}}{0.117 \times 1(\text{ml})}$$

$$= \Delta \text{OD} \times 8.547 \times \text{希釈倍率}$$

$$\text{U/mg} = \text{U/ml} \times 1/\text{C}$$

0.117 : 1mg% Purpurogallinエーテル溶液の420nmにおける吸光度

C : 溶解時の酵素濃度(c mg/ml)

(注)1Purpurogallin単位は13.5国際単位(o-dianisidineを基質とし、25℃の反応条件下)に相当する。