


● TOYOBO ENZYMES ●  
(Diagnostic Reagent Grade)

# UREASE

*from Jack bean*



## PREPARATION and SPECIFICATION

Appearance	: White amorphous powder, lyophilized
Activity	: Grade II (-201) 100U/mg-solid or more
Contaminants	: Asparaginase $\leq 2.0 \times 10^{-2}\%$ Arginase $\leq 2.0 \times 10^{-3}\%$ $NH_4^+$ $\leq 5.0 \times 10^{-4} \mu g/U$
Stabilizers	: EDTA, glutathione, succinate, BSA

## PROPERTIES

Stability	: Stable at $-20^\circ C$ for at least one year	(Fig.1)
Molecular weight	: approx. 480,000 <sup>1)</sup>	
Isoelectric point	: 5.0–5.1 <sup>1)</sup>	
Michaelis constant	: $1.05 \times 10^{-2} M$ (Urea) <sup>1)</sup>	
Structure	: 8 active sites with SH-groups per the enzyme molecule <sup>2)</sup>	
Inhibitors	: Heavy metal ions ( $Ag^+$ , $Hg^{++}$ , etc.)	
Optimum pH	: 6.0	(Fig.3)
Optimum temperature	: $60^\circ C$	(Fig.4)
pH Stability	: pH 5.5–8.5 ( $30^\circ C$ , 17hr)	(Fig.5)
Thermal stability	: below $50^\circ C$ (pH 8.0, 60min)	(Fig.6)
Effect of various chemicals	: (Table.1)	

## APPLICATIONS <sup>3)</sup>

This enzyme is useful for enzymatic determination of urea in clinical analysis.

## ASSAY

### Principle:



The disappearance of NADPH is measured at 340nm by spectrophotometry.

### Unit definition:

One unit causes the formation of two micromoles of ammonia per minute under the conditions described below.

### Method:

#### Reagents

- A. A. Urea solution : 6.0M (36g of Urea / 100ml of H<sub>2</sub>O)(Should be prepared fresh)
- B. Tris-HCl buffer, pH 8.0 : 50mM
- C.  $\alpha$ -Ketoglutarate solution : 0.25M (Dissolve 730mg of  $\alpha$ -ketoglutarate in 15 ml of H<sub>2</sub>O, adjust pH to 5.0 $\pm$ 0.1 with 5N NaOH and fill up to 20ml with H<sub>2</sub>O)(Should be prepared fresh)
- D. NADPH solution : 15mM [Dissolve 136mg of NADPH · Na<sub>4</sub> · 4H<sub>2</sub>O / 10ml of H<sub>2</sub>O] (Should be prepared fresh)
- E. Working solution (Prepare before use and store on ice)
- |       |                                  |     |
|-------|----------------------------------|-----|
| 69 ml | Tris-HCl buffer                  | (B) |
| 0.3ml | $\alpha$ -Ketoglutarate solution | (C) |
| 1.8ml | NADPH solution                   | (D) |
| 0.9ml | H <sub>2</sub> O                 |     |
- F. GIDH (glutamate dehydrogenase) solution : ca.1,000U/ml [Toyobo GradeII ,GTD-209 (Tris-HCl buffer solution, free from ammonia)]
- G. Enzyme diluent : 10mM K-phosphate buffer containing 20mM EDTA and 0.2% BSA, pH 7.0

#### Procedure

- Prepare the following reaction mixture in a cuvette (d=1.0cm) and equilibrate at 37°C for about 5 minutes.
 

2.40ml	Working solution	(E)
0.05ml	GIDH solution	(F)
0.35ml	H <sub>2</sub> O	
0.10ml	Enzyme solution*	

Concentration in assay mixture	
Tris-HCl buffer	38 mM
Urea	200 mM
$\alpha$ -Ketoglutarate	0.83mM
NADPH	0.30mM
EDTA	0.67mM
GIDH	ca.17 U/ml

- Add 0.10ml of urea solution (A) and mix by gentle inversion.
- Record the decrease in optical density at 340nm against water for 3 to 4 minutes in a spectrophotometer thermostated at 37°C, and calculate the  $\Delta OD$  per minute from the initial linear portion of the curve ( $\Delta OD$  test).  
At the same time, measure the blank rate ( $\Delta OD$  blank) by using the same method as the test except that the enzyme diluent is added instead of the enzyme solution.

\* Dissolve the enzyme preparation in ice-cold enzyme diluent (G) and dilute to 0.07–0.25U/ml with the same buffer and store on ice.

#### Calculation

Activity can be calculated by using the following formula :

$$\text{Volume activity (U/ml)} = \frac{\Delta OD/\text{min} (\Delta OD \text{ test} - \Delta OD \text{ blank}) \times V_t \times df}{6.22 \times 2 \times 1.0 \times V_s} = \Delta OD/\text{min} \times 2.41 \times df$$

$$\text{Weight activity (U/mg)} = (\text{U/ml}) \times 1/C$$

V<sub>t</sub> : Total volume (3.0ml)

V<sub>s</sub> : Sample volume (0.10ml)

6.22 : Millimolar extinction coefficient of NADPH at 340nm (cm<sup>2</sup>/micromole)

2 : Factor based on the fact that hydrolysis of one mole of urea is equivalent to oxidation of two moles of NADPH

1.0 : Light path length (cm)

df : Dilution factor

C : Enzyme concentration in dissolution (c mg/ml)


## REFERENCES

- 1) J.B.Sumner and D.B.Hand; *J.Am.Chem.Soc.*, **31**, 1255 (1929)
- 2) G.Gorin and C.C.Chin; *Biochim.Biophys.Acta.*, **99**, 418 (1965)
- 3) *Rinsho Kagaku Bunseki(Japanese)*, II ,p18(M.Kitamura,M.Saito and M.Niwa,ed.)Tokyo Kagaku Dojin,Tokyo (1969)
- 4) H.G.Schlegel and H.Kaltwasser; *Methods of Enzymatic Analysis*, Vol.2,p1081 (H.U.Bergmeyer,ed.), Verlag Chemie Weiheim, Academic Press, New York-London (1974)


**Table 1 Effect of Various Chemicals on Urease**

[The enzyme dissolved in 20mM phosphate buffer, pH 7.0 was incubated with each chemical at 30°C for 1hr.]

Chemical	Concn.(mM)	Residual activity(%)	Chemical	Concn.(mM)	Residual activity(%)
None	—	100	MnCl <sub>2</sub>	1.0	66
NaCl	10	96	MgCl <sub>2</sub>	1.0	97
Na <sub>2</sub> SO <sub>4</sub>	10	104	CaCl <sub>2</sub>	1.0	105
CH <sub>3</sub> COONa	10	108	ZnCl <sub>2</sub>	1.0	104
Na <sub>2</sub> HPO <sub>4</sub>	10	100	FeSO <sub>4</sub>	1.0	94
Citrate-Na <sub>2</sub>	10	100	CuSO <sub>4</sub>	1.0	99
Na <sub>2</sub> CO <sub>3</sub>	10	100	Ag <sub>2</sub> SO <sub>4</sub>	0.1	9
Na <sub>3</sub> BO <sub>4</sub>	10	104	HgCl <sub>2</sub>	0.1	8
Na <sub>3</sub> S <sub>2</sub> O <sub>4</sub>	10	108			




**Fig. 1. Stability (Powder form)**  
[kept under dry conditions]




**Fig. 3. pH-Activity**

[30°C in 10mM buffer solution: pH3.0-9.0]  
[Veronal-CH<sub>3</sub>COONa-HCl; pH9.0-11.0,  
glycine-NaOH.]




**Fig. 5. pH-Stability**


[30°C ,17hr-treatment with 10mM buffer  
solution: pH 3.0-9.0, Veronal-CH<sub>3</sub>COONa-HCl;  
[pH9.0-11.0: glycine-NaOH]



**Fig. 2. Stability (Powder form)**  
[kept under dry conditions]




**Fig. 4. Temperature activity**  
[in 20mM phosphate buffer,pH7.0]



**Fig. 6. Thermal stability**  
[60min-treatment with 20mM phosphate  
buffer, pH8.0.]

## 活性測定法 (Japanese)

### 1.原理



NADPHの減少量を340nmにおける吸光度の変化で測定する。

### 2.定義

下記条件下で1分間に2マイクロモルのアンモニアを生成する(1マイクロモルの尿素を加水分解する)酵素量を1単位(U)とする。

### 3.試薬

- 6.0M 尿素水溶液(36gの尿素を蒸留水に溶解し100mlとする)(用時調製)
- 50mM Tris-HCl緩衝液, pH8.0
- 0.25M  $\alpha$ -ケトグルタル酸水溶液(730mgの $\alpha$ -ケトグルタル酸を蒸留水約15mlに溶解し, 5N NaOHでpH5.0 $\pm$ 0.1に調製後, 蒸留水にて20mlとする)(用時調製)
- 15mM NADPH水溶液 [136mgのNADPH $\cdot$ Na $_2$  $\cdot$ 4H $_2$ Oを蒸留水10mlに溶解する] (用時調製)
- 試薬混液(使用直前に調製し, 氷冷保存する)
 

69.0ml	Tris-HCl緩衝液	(B)
0.3ml	$\alpha$ -ケトグルタル酸水溶液	(C)
1.8ml	NADPH水溶液	(D)
0.9ml	蒸留水	
- グルタミン酸脱水素酵素(GIDH)溶液 [東洋紡製, GTD-209(約1,000U/mlに蒸留水で希釈して使用する)]

酵素溶液: 酵素標品を予め氷冷した20mM EDTAと0.2% BSAを含む10mM K-リン酸緩衝液, pH7.0で溶解し,同緩衝液で0.07~0.25U/mlに希釈して氷冷保存する。

### 4.手順

- 下記反応混液をキュベット(d=1.0cm)に調製し,37°Cで約5分間予備加温する。
 

2.40ml	試薬混液	(E)
0.05ml	GIDH溶液	(F)
0.35ml	蒸留水	
0.10ml	酵素溶液	
- 基質溶液(A)0.10mlを添加し,ゆるやかに混和後,水を対照に37°Cに制御された分光光度計で340nmの吸光度変化を3~4分間記録し,その初期直線部分から1分間当りの吸光度変化を求める( $\Delta$ ODtest)。
- 盲検は反応混液①に酵素溶液の代わりに酵素希釈液(20mM EDTAを含む10mM K-リン酸緩衝液, pH7.0)を0.10ml加え,上記同様に操作を行って,1分間当りの吸光度変化を求める( $\Delta$ ODblank)。

### 5.計算式

$$U/ml = \frac{\Delta OD/min (\Delta OD \text{ test} - \Delta OD \text{ blank}) \times 3.0(ml) \times \text{希釈倍率}}{6.22 \times 2 \times 1.0 \times 0.10(ml)}$$

$$= \Delta OD/min \times 2.41 \times \text{希釈倍率}$$

$$U/mg = U/ml \times 1/C$$

6.22 : NADPHのミリモル分子吸光係数 (cm<sup>2</sup>/micromole)

2 : 酵素反応で1分子の尿素の加水分解に由来するNADPHの酸化は2分子である事による係数

1.0 : 光路長(cm)

C : 溶解時の酵素濃度(c mg/ml)