


● TOYOBO SUBSTRATES ●
(Diagnostic Reagent Grade)

Gal-G2- α -CNP

for Amylase substrates



SPECIFICATION

Appearance	: Slightly yellowish amorphous powder
Solubility	: Clear solution in water
Content (HPLC)	: $\geq 90\%$

PROPERTIES

Stability	: Stable at -20°C for at least one year
Molecular weight	: 659.98

APPLICATIONS

This substrate makes it possible to measure amylase activity in samples without auxiliary enzymes.