

anti-Procalcitonin monoclonal antibodies


Specifications

Specificity : human procalcitonin (PCT)
 Source : mouse
 Purification : affinity purification
 Presentation : PBS, pH7.4, 0.09% sodium azide
 Storage : -20°C (4°C for short term)

Applications

Lateral flow, ELISA, and others

Lateral flow assay using PCT-103 as the capture antibody and PCT-104 as the detection antibody



Ordering Information

Product	Catalogue No.	mAb Isotype	Usage
anti-Procalcitonin	PCT-103	IgG1 - kappa	Capture
anti-Procalcitonin	PCT-104	IgG2b - kappa	Detection