

anti-PSA monoclonal antibodies


Specifications

Specificity : human prostate specific antigen (PSA)
 Source : mouse
 Purification : affinity purification
 Presentation : PBS, pH7.4, 0.09% sodium azide
 Storage : -20°C (4°C for short term)

Applications

Lateral flow, ELISA, and others

Lateral flow assay using PSA-101 as the capture antibody and PSA-103 as the detection antibody



Ordering Information

Product	Catalogue No.	mAb Isotype	Usage
anti-PSA	PSA-101	IgG1 - kappa	Capture
anti-PSA	PSA-103	IgG1 - kappa	Detection